

# CARE PLUS

“High Dependency” Open Care Incubator



CE 1023



# CARE PLUS

## “High Dependency” Open Care Incubator

### The unique Care Plus Servo model comprises:

Three modes of warming:

- Servo Skin Control
- Air Control
- Manual Mode.

In manual mode the user can select the heater output and the time in minutes.

### The Standard features includes

- Large Baby bed with drop down acrylic side panels.
- CE marked ceramic heater for uniform warming of Baby bed ensures no Hot or Cold zones.
- Two LED based observation lamps.
- One monitor shelve to accommodate necessary monitoring devices.
- One stand provided to hold syringe pumps
- Baby bed can be tilted to  $\pm 12^\circ$  manually.
- Unit mounted on heavy duty castors. Two with locking device.
- X ray compartment provided below baby tray for trouble free X ray procedures.
- Three close drawers to store necessary baby related material.
- Rotatable Heater Box for ease in X ray procedures.



Life changes so do concepts, ideas hence for constant up gradation we reserve the right to change the specifications & features without prior notice to ensure harmony between man & the machine.



Innovated by

## Shreeyash Electro Medicals

Sr. No. 49, Plot No. 22, Shri Sai Industrial Estate,  
Gujarwadi Phata, Katraj, Pune 411046, INDIA

Mobile: +91 98223 75996 / +91 94223 21885

[www.shreeyashindia.com](http://www.shreeyashindia.com)



ISO 13485 : 2008

CE 1023