



CARE PLUS

“High Dependency” Open Care Incubator



CARE PLUS

“High Dependency” Open Care Incubator

The unique Care Plus Servo model comprises:

Three modes of warming:

- Servo Skin Control
- Air Control
- Manual Mode.
In manual mode the user can select the heater output and the time in minutes.
- A separate LCD display provided to inform the user regarding various messages & errors.

The Standard features includes

- Large Baby bed with drop down acrylic side panels.
- CE marked ceramic heater for uniform warming of Baby bed ensures no Hot or Cold zones.
- Two LED based observation lamps.
- One sliding monitor shelve provided to accommodate necessary monitoring devices.
- One large stand provided to hold syringe pumps
- Smooth and noiseless tray tilt mechanism. Baby bed can be tilted upto $\pm 12^\circ$ manually.
- Unit mounted on heavy duty castors. Two with locking device.
- X ray compartment provided below baby tray for trouble free X ray procedures.
- Three close drawers to store necessary baby related material.
- Rotatable Heater Box for ease in X ray procedures.
- One I.V. rod is a standard accessory.



Life changes so do concepts, ideas hence for constant up gradation we reserve the right to change the specifications & features without prior notice to ensure harmony between man & the machine.



Innovated by

Shreeyash Electro Medicals

Sr. No. 49, Plot No. 22, Shri Sai Industrial Estate,
Gujarwadi Phata, Katraj, Pune 411046, INDIA
Mobile: +91 98223 75996 / +91 98224 76056
Factory: +91 89561 44211 / +91 89561 44212
dragonfly.sudhir@gmail.com | www.shreeyashelctromedicals.com

ISO 13485 certified
CDSCO registered