

INCUBATOR





Servo Deluxe Incubator

- The extra large baby bed that will softly comfort and effectively treat even a full term infant
- The withdrawable baby tray with auto lock facility, offers greater access to infant for delicate procedures
- The externally adjusted head up/down position ensures undisturbed and calm baby chamber environment.
- The servo proportional incubator offers "Exclusive, Innovative, Never before in the medical world", incorporation of 3 level safety circuit for absolute and maximum safety of the infant
- The precisely manufactured servo proportional controller maintains pin point accuracy of environment and baby temperature.
- The L.E.D. bargraph indicates the precise level of heater output at any given time.
- The detachable sub assemblies and auto fumigation mode ensures clean sterile ambience in the infant canopy.
- The aerodynamic styling of incubator saves space and offers maximum dexterity, to the use
- A two way air distributor system ensures a uniformity of temperature in baby chamber and avoids radiation heat losses.
- A separate humidity booster with control and preassembled hygrometer ensures easy control and indication of relative humidity in the baby compartment.



Transport Incubator

Comprises:

- Double wall canopy to avoid radiative heat losses.
- Forced air circulation for uniform air temp in canopy.
- Micro controller based Temp. Indicator & controller along with 4 Level Safety System.
- Baby tray provided with head low, high facility.
- Front flap with two nursing ports.
- Baby holding straps provided for infant fixation during transit.

The Thoughtful Design of the Trolley Facilitates:

- Horizontal fixation of the entire system in small ambulance behind driver's seat.
- Fits within the smallest elevator (size 38 (L) x 39 (B))
- Electrically operated neonatal slow suction (optional at extra cost)
- 2 Handles on each side of trolley for easy transportation

- In line water Nanometer to indicate gas pressure delivered to Baby Via Ambu bag (optional at extra cost)
- Trolley incorporated with necessary power pack along with 2 maintenance free batteries.
- Incubator operates on battery for 3 hours (180 minutes) in fully charged battery condition.



Fluid Warming Cabinet

- Comprises forced air warming system for temperature in the cabin.
- Micro controller based servo proportionate system incorporated to maintain pin point accuracy in the chamber.
- Special platinum PT-100 Sensor with steel capsule.
- Digital display provided to indicate set and actual temperature
- Accommodates 12, 500 ML fluid bottles.
- Front door provided with glass for easy viewing.
- Durable CE marked ceramic heater which has incredible long life.
- Adjustable shelf facility.
- Illumination lamp provided in the cabinet.

Life changes so do concepts, ideas hence for constant up gradation we reserve the right to change the specifications & features without prior notice to ensure harmony between man & the machine.



Innovated by

Shreeyash Electro Medicals

Sr. No. 49, Plot No. 22, Shri Sai Industrial Estate,

Gujarwadi Phata, Katraj, Pune 411046, INDIA

Mobile: +91 98223 75996 / +91 98224 76056

Factory: +91 89561 44211 / +91 89561 44212

dragonfly.sudhir@gmail.com | www.shreeyashelectromedicals.com

ISO 13485 certified
CDSCO registered